



Bone Marrow Processing

Advancing Therapeutic Apheresis and Cell Collections to the Next Level of Patient Care

Bringing ease of use to bone marrow processing procedures

The Spectra Optia system helps simplify the bone marrow processing (BMP) procedure through on-screen guidance and a single-bag approach. With a minimal amount of interaction needed throughout the procedure, you can perform other tasks while the system processes the bone marrow in under 90 minutes.

Procedure and System Highlights


Minimal required interactions	<ul style="list-style-type: none">Streamlines procedure with a single-bag approachRequires less operator supervision, as the system continuously monitors and adjusts the interfaceMinimizes memorization and simplifies training
Consistent results	<ul style="list-style-type: none">Automated interface management (AIM) system is designed to produce consistent results through interface stabilityPerforms red blood cell (RBC) reduction greater than 97%*Provides volume reduction greater than 90%, producing a final product volume under 200 mL*
Intuitive graphical user interface (GUI)	Streamlines your procedure management with on-screen instructions and touch-screen data entry
Procedural flexibility	Provides you with control over parameters to meet your target processing outcomes
Large, durable wheels on pivoting casters	Makes maneuvering the system throughout your facility easy

*Validation report on file at Terumo Blood and Cell Technologies.



How It Works

The GUI guides you through each step of the procedure. The examples below demonstrate the ease of monitoring and controlling the settings throughout a procedure.



BMP connection screen


The system guides you through connecting the single BMP bag to the inlet and return line on the IDL Set.

Collection status screen

Continuously compares how the system is collecting cells with the target collection preference

Gives you the ability to adjust your collection preference

Displays an image of the collect port throughout your collection procedure



	Target	Current
BM Processed (mL)	6000	6000
BM Cycles	4.0	4.0
Collection Bag (mL)	100	100
Plasma Bag (mL)	100	100

Run targets attained.

16:59
19-10-2012


Rinseback

BMP

Run target screen

Once your run targets are attained, the system sounds a tone and displays the target and current values. Yellow frames around a current value make it immediately apparent which run targets were attained.

BMP Accessory Set*



*Must be used in conjunction with the Spectra Optia® IDL Set.

Compact packaging	Minimizes the space required for storage
Color-coded components	Simplify setup and operation

Working With You

Every interaction we have with you is important. By fostering open and ongoing relationships, we bring more value to you and the patients we are all focused on serving.

Even after the technology is in place, we continue to serve you with:

- Education and training
- Technical support
- Clinical and scientific support
- Customer support
- User groups and professional networks

Product and procedure availability varies according to regulatory clearance or approval in each country.



Terumo Blood and Cell Technologies is a medical technology company. Our products, software, and services enable customers to collect and prepare blood and cells to help treat challenging diseases and conditions. Our employees around the world believe in the potential of blood and cells to do even more for patients than they do today. [TerumoBCT.com](https://www.terumobct.com)

Terumo BCT, Inc. Lakewood, CO, USA +1.303.231.4357	Terumo BCT Europe N.V. Zaventem, Belgium +32.2.715.0590	Terumo BCT Asia Pte. Ltd. Singapore +65.6715.3778	Terumo BCT Latin America S.A. Buenos Aires, Argentina +54.11.5530.5200	Terumo BCT Japan, Inc. Tokyo, Japan +81.3.6743.7890
----------------------------------------------------------	---------------------------------------------------------------	---------------------------------------------------------	------------------------------------------------------------------------------	-----------------------------------------------------------